STORZ

KARL STORZ—ENDOSKOPE

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C-MAC® Videolaryngoskop Serie 8403 xx, Anschlusskabel 8403 X

INSTRUCTION MANUAL

C-MAC® Video laryngoscope series 8403 xx, Connection cable 8403 X

MANUAL DE INSTRUCCIONES

Videolaringoscopio C-MAC® serie 8403 xx, cable de conexión 8403 X

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It is recommended to check the suitability of the products for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper use, cleaning and sterilization of the C-MAC® video laryngoscope and the monitor 8403 ZX. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

WARNING: KARL STORZ instruments are provided unsterilized, and must be cleaned and disinfected or sterilized prior to the initial use and before each subsequent use.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed physician.

NOTE: All users must have been instructed in the operation and use of the C-MAC® video laryngoscopes including the C-MAC® Monitor 8403 ZX or C-MAC® PM 8403 XD, otherwise this monitor may not be used.

NOTE: The C-MAC® video laryngoscopes including the C-MAC® Monitor 8403 ZX or C-MAC® PM 8403 XD must be used in accordance with the airway algorithm in the respective country of application or the official specifications for airway management.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be deduced from the information, figures and descriptions in this manual.

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1	Controls	4
2	Symbol description.....	6
2.1	Pictograms	6
2.2	Symbols for safety instructions	7
3	Intended use	8
4	User qualification.....	8
4.1	Training in the operation and function of the device	9
4.2	Patient population.....	9
4.3	User profile of physician and assistants	9
5	Safety Instructions	10
5.1	Operation.....	13
5.1.1	Connection to the C-MAC® monitor.....	14
5.2	Connection to the C-MAC® PM	15
5.2.1	Switching the OTI display on	16
5.2.2	Switching the OTI display off.....	16
5.3	Connection to the C-HUB® II	17
5.4	Multifunctional button of the video laryngoscope	18
5.4.1	Color coding meanings for the multifunctional button 'BlueButton'	18
5.4.2	Single image capture	18
5.4.3	Video recording	18
5.4.4	Inserting an oxygen or suction catheter (only with 8403 AX/BX/HX)	19
5.5	Disassembly	19
6	Reprocessing (worldwide except US and Canada)	20
6.1	Overview table	20
6.2	General Warnings	21
6.3	Accessories	22
6.4	Manual wipe-down disinfection.....	22
6.5	Preparation for cleaning and disinfection	23
6.6	Manual cleaning.....	23
6.7	Manual disinfection	23
6.8	Machine cleaning and disinfection	23

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


Controls

- ① C-MAC® video laryngoscope blade
- ② Socket for Connection cable ④
- ③ Multifunctional button (BlueButton)
- ④ Connection cable
- ⑤ Socket to the monitor
- ⑥ Socket to the laryngoscope
- ⑦ C-MAC® pocket monitor
- ⑧ TFT display
- ⑨ Coupling element
- ⑩ C-MAC® monitor
- ⑪ Connection 2 for video endoscope
- ⑫ Connection 1 for video endoscope
- ⑬ Mains connection
- ⑭ Reset button
- ⑮ Holder for SD memory card
- ⑯ Status bar
- ⑰ Pilot lamp
- ⑱ On/Off switch
- ⑲ Pilot lamp charging (orange)
- ⑳ Softkeys
- ㉑ Function symbols
- ㉒ USB port
- ㉓ HDMI port

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Avoidance of environmental pollution through the use of electronic equipment (China RoHS)

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WEEE symbol: Labeling of electrical devices

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Labeling of batteries/accumulators

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UL approval mark, valid in the USA and Canada

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Degree of protection in accordance with International Protection Codes: Protection against contact with tools or similar and small foreign bodies, protection against vertically dripping water

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
Rated ambient temperature/
transducer temperature class

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2.2 Symbols for safety instructions

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

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3 Intended use

C-MAC® Video Laryngoscopes:

The C-MAC® Video Laryngoscopes 8403 AX/BX/HX/AXC/BXC/DXC/EXC/GXC/KXC/MXC/HXP is used in combination with Connection cable 8403 X and Monitor 8403 ZX or the C-MAC® PM 8403 XD for endotracheal intubation and the inspection of the oropharynx.

C-MAC® Monitor:

The C-MAC® Monitor is used together with a CMOS videoscope to play back images during an endoscopic intervention.

Contraindications: No contraindications relating directly to the product are currently known. Use of the C-MAC® video laryngoscope is contraindicated if, in the opinion of the responsible physician, the health of the patient is endangered through its use, e.g., due to the general condition of the patient, or if laryngoscopy itself is contraindicated.

4 User qualification

The C-MAC® video laryngoscope may only be used by persons with an appropriate medical qualification and who are acquainted with the laryngoscopic technique. The information given in these instructions only serves to instruct in the correct handling, cleaning, and sterilization of the C-MAC® video laryngoscope. It is not intended as an introduction to laryngoscopy.

The C-MAC® video laryngoscopes must be used in accordance with the airway algorithm in the respective country of application or the official specifications for airway management.

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4.1 Training in the operation and function of the device

Your local representative or responsible KARL STORZ member of staff is available to provide training and information on further training alternatives.

4.2 Patient population

Men and women as well as adolescents.

State of health:

suitable for laryngoscopic interventions, in the physician's judgement.

4.3 User profile of physician and assistants

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in one of the languages used on the device and in the instruction manual.
- Thorough training in the operation of the instruments.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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5 Safety Instructions

Risks of infection

- Warning: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

- Warning: Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

Risks due to electric energy

- Warning: When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. This is particularly important if a type CF applied part endoscope is used, in which case a type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.

Risks due to damaged parts

- Warning: Check the C-MAC® video laryngoscopes prior to each use to ensure that they function correctly and have been properly cleaned, disinfected and/or sterilized.
- Warning: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners,

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burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

- Warning: If damage is apparent, the C-MAC® video laryngoscope should no longer be used, if necessary it should be sent off to be repaired.
- Warning: C-MAC® blades contain a nickel compound and can trigger a nickel allergy in allergy-prone people.

Risks due to combination of system components

- Warning: Combinations of medical devices are only then assured to be safe if
 - they are identified as such in the respective instruction manual or
 - the intended purpose and the interface specifications of the devices used in combination permit this.
 - Warning: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. All relevant instruction manuals must be read carefully and the instructions given must always be followed precisely. Check the functioning of the products used in combination.
 - Warning: Additional equipment connected to medical electrical equipment must be able to prove their compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively).
- Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned

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standard requirements. If in doubt, please consult your local specialist dealer or the technical service department (Standard/directive references: IEC 60601-1+A1+A2:1995: 6.8.2.c, 19.2.b, 19.2.c, IEC 60601-1:2005: 7.9.2.5, 8.1, 16.2.d, MDD 93/42/EEC: Annex I clause 13.6.c).

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Damage to the instrument

- Caution: The C-MAC® video laryngoscopes must not be cleaned in an ultrasound bath.
- Caution: The C-MAC® video laryngoscopes must not be steam-sterilized.

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5.1 Operation

The C-MAC® video laryngoscopes 8403 xx can be used together with the following monitors:

- C-MAC® monitor 8403 ZX (Software version 704v300 or higher)
- C-MAC® PM Pocket Monitor 8403 XD
- with additional monitors via the camera control unit C-HUB® II 20 2903 01 (Software version: 721v105 or higher)

 **WARNING:** Before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the endoscope provides a live image (rather than a stored one) and has the correct image orientation

 **NOTE:** You can remove the connection cable from the video laryngoscope when in operation and insert it in another laryngoscope.

 **WARNING:** It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other units at any time.

 **WARNING:** Risk of injury: Optical radiation can lead to eye injuries. Never look into the light exit at the tip of the video laryngoscope.

 **NOTE:** To prevent the telescope from fogging up, the C-MAC® video laryngoscope should not be cooled below 18 °C and should be switched on in good time to allow the optical system to warm up. In extreme situations, the antifog solution "ULTRA STOP" with article number 15006 B can also be used.

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5.1.1 Connection to the C-MAC® monitor

NOTE: The video laryngoscopes 8403 xx can only be connected to the C-MAC® monitor 8403 ZX. Older monitor types 8402/8401 ZX are not compatible.

NOTE: The instruction manual of the C-MAC® monitor must be observed (Art. no. 96076008D)!

1. Connect the Connection cable ⑥ to the socket ② of the video laryngoscope. The orientation pin on the plug and socket facilitates connection.

NOTE: Please make sure that the connection cable is fully inserted. No gap should be visible.

2. Connect the Connection cable ⑤ coming from the video laryngoscope to the connection socket ⑫ (primary connection) or the connection socket ⑪ of the C-MAC® monitor.

NOTE: The monitor 8403 ZX enables two CMOS video units to be connected with Lemo connectors. All CMOS videoscopes from KARL STORZ can be connected to both connections. Video units from other manufacturers cannot be connected.

NOTE: If two cameras are connected, the user must select the camera for which the images are to be displayed in the menu (refer to the instruction manual for the C-MAC® monitor 8403 ZX, Art. no. 96076008D).

3. Switch on the monitor by pressing the power switch ⑮. The green pilot lamp ⑮ on the monitor and the blue multifunctional button ③ on the video laryngoscope light up. The devices are now ready for use.

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5.2 Connection to the C-MAC® PM

NOTE: The video laryngoscopes 8403 xx can only be connected to the C-MAC® PM 8403 XD. The 8401 XD model is not compatible.

NOTE: The instruction manual of the C-MAC® PM 8403 XD must be observed (Art. no. 96076020D)!

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1. Connect the C-MAC® PM with the socket ② of the video laryngoscope blade. The orientation pin on the plug and socket facilitates connection.

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NOTE: The monitor can also be opened up when connecting both components.

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2. The devices are now ready for use. The multifunction button ③ of the video laryngoscope lights up blue and indicates that still image capture and the video function are ready for use.

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NOTE: When inserting the C-MAC® PM into the video laryngoscope always ensure that the lugs of the plug connection are opposite each other.

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5.2.1 Switching the OTI display on

To connect the C-MAC® PM to the video laryngoscope, the following steps must be performed:

1. Flip the OTI display ⑧ upward (photo B, C).
2. Flip the OTI display to the left (photo D).
3. Set the OTI display so that you can see the image optimally (photo E).

- NOTE:** The image focus is adjusted automatically. Manual adjustment is not necessary.
- NOTE:** Depending on the composition of the image, there may be a slight image reduction in the first 5-8 seconds.

The C-MAC® video laryngoscope with C-MAC® PM is ready for operation.

5.2.2 Switching the OTI display off

1. Put the OTI display ⑧ into the vertical position (photo D).
2. Flip the OTI display to the right (photo C).
3. Flip the OTI display downward (photo B).
4. The C-MAC® PM is now turned off.

- NOTE:** You can remove the C-MAC® PM from the C-MAC® video laryngoscope when in operation and insert it in another C-MAC® video laryngoscope without having the C-MAC® PM turned off.

WARNING: After 10 minutes, the display of the C-MAC® PM turns off automatically and **must be turned back on again by closing it and reopening it** (switching it off and on again).

- NOTE:** It is not necessary to perform a white balance.

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5.3 Connection to the C-HUB® II


NOTE: The video laryngoscopes 8403 xx can only be connected to the C-HUB® II 20 2903 01 (software version 721v105 or higher). Older software versions are not compatible..

NOTE: The instruction manual of the C-HUB® II must be observed (Art. no. 96206234D)!

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1. Connect the Connection cable ⑥ with the socket ② of the video laryngoscope. The orientation pin on the plug and socket facilitates connection.

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2. Connect the power supply unit to the C-HUB® II.
3. Connect the Connection cable ⑤ coming from the video laryngoscope to the video input on the front of the C-HUB® II.

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4. Connect the monitor to the S-Video or HDMI output on the rear of the C-HUB® II.

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5.4 Multifunctional button of the video laryngoscope

Both single image capture and video recording can be activated using the multifunctional button ③ on the video laryngoscope.

Documentation with the multifunctional button of the video laryngoscope can be performed on the C-MAC® monitor 8403 ZX and C-MAC® PM 8403 XD.

Approx. 5 seconds after each time the C-MAC® PM is switched on, a blue bar is displayed at the top of the screen for a few seconds. The remaining memory capacity in minutes is shown on the left-hand side, while the name of the directory in which the saved data are stored is shown on the right.

5.4.1 Color coding meanings for the multifunctional button "BlueButton"

As soon as the multi-functional button on the C-MAC® video laryngoscope blade lights up blue, a blue information bar will appear on the TFT monitor for a few seconds.

It is now possible to record images and videos.

The colors of the multi-functional button have the following meanings:

Blue: permanently illuminated The recording function is enabled

Green: temporarily illuminated when recording an image, or flashes during video recording.

5.4.2 Single image capture

Single image capture is activated when the multifunctional button ③ is pressed down once briefly. The current image on the screen is saved on the SD memory card inserted into the connected monitor 8403 ZX.

During single image capture, the multifunctional button ③ briefly lights up green and then switches back to blue.

5.4.3 Video recording

Video recording is started by pressing and holding the multifunctional button ③ (for approx. 2 seconds) and is stopped by briefly pressing the multifunctional button once more. The video

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recording is saved on the SD memory card inserted into the connected monitor 8403 ZX. During video recording, the multifunctional button ③ on the video laryngoscope emits a pulse-like green light which switches back to blue when pressed again briefly at the end of recording.

NOTE: The video stream can only be viewed on a PC if an MPEG 4-codec is installed on the PC.

NOTE: If memory capacity is less than 10 minutes, the bar is displayed in orangered.

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5.4.4 Inserting an oxygen or suction catheter (only with 8403 AX/BX/HX)

When using a blade with an oxygen or suction catheter guide, the probe or catheter is advanced so far beyond the guide until visibility in the image is guaranteed.

When using the blade 8403 AX, the 14 Fr. (green) and the 16 Fr. (orange) catheters can be used.

When using the blade 8403 HX and 8403 BX, the 16 Fr. (orange) and the 18 Fr. (red) catheters can be used.

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5.5 Disassembly

Disconnect the connecting cord from the C-MAC®. Disconnect the connecting cord or the C-MAC® PM from the videolaryngoscope.

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**6 Reprocessing (worldwide
except US and Canada)**

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- This reprocessing step is part of the validated process
 - This reprocessing step may be carried out additionally in order to aid cleaning
 - X** On the grounds of potential damage, this reprocessing step must NEVER be performed
- No symbol:** This reprocessing step has neither been validated in respect of efficacy nor in respect of material compatibility by KARL STORZ.

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Reprocessing (worldwide except US and Canada)

6.2 General Warnings
WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

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NOTE: As a general rule, all system components of the 8403 xxx C-MAC® video laryngoscope and the C-MAC® connecting cable 8403 X are resistant to manual and machine cleaning and disinfection. If the protective cap is not being used, then manual rinsing and drying of the contacts may be necessary in order to ensure that functioning is not impaired. For this reason, the protective cap 8403 YZ is also offered. The user is responsible for deciding whether or not this needs to be used, depending on the individual process management (drying).

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

NOTE: For information on reprocessing of the C-MAC® Monitor 8403 ZX, please refer to the corresponding instruction manual 96076008D.

NOTE: For information on reprocessing of the C-MAC® PM, please refer to the corresponding instruction manual 96076020D.

6.3 Accessories

Accessories required for carrying out reprocessing:

- Brushes: 27652

6.4 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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i **NOTE:** Validation was performed with the Tristel Trio Wipes System from Tristel Solutions Limited.

6.5 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

6.6 Manual cleaning

The medical devices must be completely immersed in a cleaning solution. Through opening and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. After the necessary exposure time, clean the instrument mechanically with the aid of a brush or sponge. Finally, it must be rinsed with cold water to ensure neutralization.

6.7 Manual disinfection

The medical devices must be completely immersed in a disinfectant solution. Through opening and filling the lumina in a targeted manner, it must be ensured that even surfaces with restricted access are covered and no air bubbles are present. Following the exposure time, the medical devices must be rinsed several times with completely demineralized water or microbiologically pure, sterile water in order to remove all chemical residues. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with sterile compressed air. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

6.8 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

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Automated cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off afterwards by hand.

i NOTE: Drying was validated at 110°C.

6.9 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

i NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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6.11 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and following consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Hydrogen peroxide (H₂O₂) sterilization – ASP STERRAD®

 **WARNING:** Medical devices must be dismantled for sterilization.

 **WARNING:** Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

 **CAUTION:** For more detailed information, please consult the user handbook of the respective device.

 **NOTE:** Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

 **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.

The following STERRAD® sterilization processes have been validated and approved by KARL STORZ for this medical device:

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Reprocessing
(worldwide except US
and Canada)

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7 Reprocessing (for US and Canada customers only)

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- This reprocessing step is part of the validated process
- This reprocessing step may be carried out additionally in order to aid cleaning
- X** On the grounds of potential damage, this reprocessing step must NEVER be performed
No symbol: This reprocessing step has neither been validated in respect of efficacy nor in respect of material compatibility by KARL STORZ.

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7.2 General Warnings

 **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction.

 **WARNING:** Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.

 **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

 **CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage.

 **CAUTION:** Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing.

 **CAUTION:** National laws and regulations must be observed.

 **NOTE:** For information on reprocessing of the C-MAC® Monitor 8403 ZX, please refer to the corresponding instruction manual 96076008D.

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Reprocessing

for US and Canada
customers only

- ① **NOTE:** For information on reprocessing of the C-MAC® PM, please refer to the corresponding instruction manual 96076020D.
- ① **NOTE:** As a general rule, all system components of the 8403 xxx C-MAC® video laryngoscope and the C-MAC® connecting cable 8403 X are resistant to manual and machine cleaning and disinfection. If the protective cap is not being used, then manual rinsing and drying of the contacts may be necessary in order to ensure that functioning is not impaired. For this reason, the protective cap 8403 YZ is also offered. The user is responsible for deciding whether or not this needs to be used, depending on the individual process management (drying).

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- Brushes: 27652

7.4 Water Quality Requirements Utility Water:

Water as it comes from the tap that meets specifications as defined per AAMI TIR34:2014 Water for Reprocessing of Medical Devices. Utility water is acceptable for pre-cleaning and cleaning of the instruments.

Critical Water:

Water that is extensively treated (usually by a multistep treatment process that could include a carbon bed, softening, DI and RO, or distillation) to ensure that the microorganisms and the inorganic and organic material are removed from the water; a final submicron filtration could also be part of the treatment process. Critical water is recommended for rinsing after high-level disinfection.

7.5 Pre-cleaning

1. Immediately after a procedure, at the point of use (e.g. operating room), wipe the device with a soft, lint-free disposable cloth moistened with water or diluted mild/neutral pH enzymatic cleaning solution to remove all heavy soiling.
2. Place the device in a tray or container and transport them to the decontamination area per your institution's policy.

NOTE: If immediate cleaning is not feasible, soak the instruments in a diluted mild/neutral pH enzymatic cleaning solution. Do not soak for longer than 30 minutes.

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7.6 Manual Cleaning

1. Thoroughly rinse the device using a commercial overhead sprayer nozzle with cold water to remove all gross debris. Rinse for a minimum of 2 minutes.
2. Completely immerse the device in a diluted mild/neutral pH enzymatic cleaning solution (e.g. Enzol). Keep immersed for a minimum of 1 minute.
3. While immersed, use a soft bristle brush (P/N: 27652/3) to brush all the surfaces of the device. Brush for a minimum of 1 minute.
4. Use a low-lint wipe to wipe the entire length of the connecting cable.
5. After the immersion and brushing period, remove the device from the cleaning solution and rinse using a commercial overhead sprayer nozzle with cold water. Rinse for a minimum of 2 minutes.
6. Dry the device and the contacts with a soft, lint-free cloth or dry, oil-free compressed air (<5 psi).

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CAUTION: V-PRO sterilization may cause cosmetic changes to the camera head that do not necessarily impact the camera's functionality.

NOTE: For sterilization with gas/Steris Amsco V-Pro 1, the nationally applicable laws and regulations must be observed.

NOTE: Please use the protective cap 8403 YZ (see Fig.) for sterilization with Steris Amsco V-Pro 1.

1. Place the camera head in a FDA-cleared sterilization tray.
2. Double wrap the tray with polypropylene wrap.
3. Start the sterilizer using the appropriate cycle.

7.11 Ethylene Oxide Sterilization (EO)

1. Place the camera head in an FDA-cleared sterilization tray.
2. Double wrap the tray with polypropylene wrap.
3. Start the sterilizer using the following parameters:

Conditioning Parameters:

Temperature: 55°C (131°F)
Humidity: 70% RH
Conditioning Dwell Time: 30 minutes

Sterilization Parameters:

Sterilant: 100% Ethylene Oxide
Temperature: 55°C (131°F)
Gas Concentration: 735 ± 30 mg/l
Exposure Time: 180 minutes

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7.12 Chemical low-temperature sterilization with peracetic acid – STERIS System 1E

- ❗ **NOTE:** For sterilization with Steris System 1E, the nationally applicable laws and regulations must be observed.
- ❗ **NOTE:** Please use the protective cap 8403 YZ (see Fig.) for sterilization with Steris System 1E.

1. Place the camera in one of the following trays:
 - C1200E SS1E General Processing Tray/ Container
 - C1220E SS1E Directed Flow Processing Tray/ Container
 - C1140E SS1E Flexible Endoscope Processing Tray/Container'
 - C1160E SS1E Universal Flexible Processing (Open) Tray
2. Load the tray into the SS1E and start the system.

7.13 High-Level Disinfection

⚠ **WARNING:** High-level disinfection should ONLY be used for instruments which come into contact with intact mucous membranes.

❗ **CAUTION:** Before disinfection, the device must be thoroughly cleaned, rinsed, and dried.

❗ **CAUTION:** Hydrogen peroxide based disinfectants may cause cosmetic changes to the devices that do not impact the functionality of the devices.

❗ **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time, and service life.

This device may be chemically disinfected using CIDEX, a 2.4% glutaraldehyde solution, CIDEX OPA, a 0.55% ortho-phthalaldehyde solution, or Resert XL HLD, a 2.0% concentration of hydrogen peroxide.

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Reprocessing for US and Canada customers only

1. Prepare the disinfecting solution per manufacturer's instructions.
2. Place the device into the disinfection solution. Ensure that no air bubbles are present on the surface of the device.
3. Use the following KARL STORZ validated conditions to achieve high-level disinfection:
 - CIDEX: Immerse for 45 minutes at 25°C (77°F)
 - CIDEX OPA: Immerse for 12 minutes at 20°C (68°F)
 - Resert XL HLD: Immerse for 8 minutes at 20°C (68°F)
4. After disinfection is complete, remove the device from the disinfection solution and completely immerse in distilled/demineralized water. Keep immersed for 1 minute. Discard the water and repeat with fresh water for a total of 3 immersion rinses.

i NOTE: If using Resert XL HLD, only one immersion rinse cycle is necessary.

5. Dry the device and the contacts with a soft, lint-free, sterile cloth or clean, dry, oil-free compressed air (<5 psi).

7.14 Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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8 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

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9 Limitation of liability

As the supplier of these medical devices we only regard ourselves responsible for the safety, reliability and performance of the products if assembly, system expansions, readjustments, modifications, or repairs have been performed by people authorized by KARL STORZ, and the device has been used in accordance with its operating instructions at all times.

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10 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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11 Directive compliance

In accordance with the Medical Device Directive (MDD):

Medical device Class I.
This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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(for 8403 AX, BX, HX, AXC, BXC, DXC, EXC, GXC, KXC, HXP, MXC, X, XD, ZX)

In accordance with EN 60601-1, EN 60601-2-18, UL 2601, CSA 601.1:

- Degree of protection against electric shock:
Applied part of type BF

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13 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its service life, dispose of the device as electronic waste. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information about your local collection point. Within the scope of application of this directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

To dispose of C-MAC® laryngoscope blades for single use, no special measures are necessary. National laws and regulations must be observed.

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14 Repair program

To bridge the repair period, you will receive a unit on loan, which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
78532 TUTTLINGEN

In other countries please contact your local KARL STORZ branch or authorized dealer.

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15 Important information

For hygiene reasons or to prevent infection, telescopes, instruments and units must be cleaned, disinfected/sterilized before they are sent for repair.

We reserve the right to return contaminated instruments/units to the sender.

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16 Technical data

**Blade 8403 AX/BX/HX/AXC/BXC/DXC/EXC/
GXC/KXC/MXC/HXP**

Protection against submersion	IPX8
Camera Technology	CMOS
Camera resolution	640 × 480
Illumination	LED, white, 1 W

Operating/storage conditions for C-MAC®

Operating conditions:	0 °C – +40 °C
Storage conditions:	30 % – 70 % rel. humidity
	-20 °C – +50 °C
Optimal operating conditions:	5 % – 95 % rel. humidity
	18 °C bis 40 °C

Charging station 8403 XDL

Classification:	IEC Class II
Line frequency:	DC
Voltage:	5V
Power consumption:	2 VA

Weight:	approx. 230 g
---------	---------------

Dimensions

Wall power supply:	ca. 60 x 70 x 35 mm (without power supply adaptor)
--------------------	---

Charging station:	ca. 103 x 75 x 25 mm
-------------------	----------------------

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17 Troubleshooting

WARNING: Always disconnect the device from the mains before carrying out any maintenance work!

Error description:

Complete failure of device. Possible lighting up of status display.

Possible causes:

- Power supply failure.
- Loose connection between AC adapter plug and socket on device.
- System error

Remedy:

- Check the power supply.
- Firmly push power cord connector into socket on device.
- Reset.

Error description:

No image, TFT screen is dark.

Possible causes:

- Connection cable/C-MAC® PM not fully inserted
- Defective camera electronics.
- TFT screen is defect.

Remedy:

- Insert connection cable/C-MAC® PM as far as it will go.
- Send electronic module/C-MAC® PM or supply unit to KARL STORZ for repair.

Error description:

Cloudy picture, stripes, streaks, or similar.

Possible causes:

- Telescope of the laryngoscope is soiled.
- Contacts of connection cable/C-MAC® PM are soiled.

Remedy:

- Clean with cotton swab and alcohol solution or special cleaning paste.
- Clean the contacts of the connection cable/C-MAC® PM.

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Color distortions.

Possible causes:

- White balance not carried out correctly.
- Electronic module defect.

Remedy:

- Carry out new white balance.
- Send supply unit to KARL STORZ for repair.

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- Camera connection cable is defect.

Remedy:

- Have a new camera connection cable fitted.

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Image cannot be saved.

Possible causes:

- No memory card inserted
- Memory card full
- Memory card is not recognized

Remedy:

- Insert memory card
- Change memory card
- Insert the memory card in the switched-off C-MAC® monitor and turn the monitor on again.
- Perform reset
- Format the memory card on the PC with FAT (standard).

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Abhilfe:

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Error description:

Video stream cannot be played back on the PC monitor.

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SD memory card cannot be written on.

Possible causes:

- SD memory card is formatted incorrectly

Remedy:

- Format the SD memory card on the PC with FAT (standard).

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
Forceps designed by MAGILL, modified by BOEDEKER, length 25 cm, suitable for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 2-4

Forceps designed by MAGILL, for children, modified by BOEDEKER, length 20 cm, suitable for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 1 and 2

Bag for intubation set -C22-, ULM model, made from water resistant, hard-wearing material, washable, the bag comprises two compartments with several storage options for C-MAC® video laryngoscopes with C-MAC® pocket monitor and for conventional laryngoscopes, for use with C-MAC® pocket monitor 8403 XD, C-MAC® video laryngoscopes and conventional laryngoscopes

Protection Cap, for the C-MAC® system Interface on the C-MAC® video laryngoscope 8403 xxx and C-MAC® PM 8403 XD as well as the C-MAC® Connecting Cable 8403 X, to protect plug contact during reprocessing, cap is reusable, for use with: C-MAC® Connecting Cable 8403 X, C-MAC® PM 8403 XD, C-MAC® Video Laryngoscopes 8403 xxx, Electronic Modules 8401 X / 8402 X and Flexible Intubation Video Endoscope (FIVE)

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseas.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzucuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skårholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com
KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Høfte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraska 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ Endoscopy India Private Limited
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com